

=====

Sequence Listing could not be accepted due to errors.

See attached Validation Report.

If you need help call the Patent Electronic Business Center at (866) 217-9197 (toll free).

Reviewer: Anne Corrigan

Timestamp: [year=2008; month=8; day=27; hr=16; min=4; sec=57; ms=237; ]

=====

\*\*\*\*\*

Reviewer Comments:

SEQUENCE LISTING

<110> TransMIT Gesellschaft für Technologietransfer mbH

<120> Animal Species-specific and quantitative Detection of Central  
Nervous System Tissue in Meat  
and Meat Products

<130> TM006/Bü/Sc/Ab

<140> 10584270

<141> 2008-08-23

<150> DE 103 61 489.3

<151> 2003-12-23

—

Please remove the foreign accent mark in the above <110> response:  
foreign accent marks are non-ASCII characters, which cannot be  
processed.

The above <120> response exceeds the Sequence Rules' required 72-  
character line limit. Please adjust the line.

\*\*\*\*\*

COPY

Application No: 10584270

Version No: 1.0

Input Set:

Output Set:

Started: 2008-08-23 06:00:55.677  
Finished: 2008-08-23 06:00:57.468  
Elapsed: 0 hr(s) 0 min(s) 1 sec(s) 791 ms  
Total Warnings: 6  
Total Errors: 6  
No. of SeqIDs Defined: 6  
Actual SeqID Count: 6

Error code	Error Description
W 213	Artificial or Unknown found in <213> in SEQ ID (1)
E 257	Invalid sequence data feature in <221> in SEQ ID (1)
W 213	Artificial or Unknown found in <213> in SEQ ID (2)
E 257	Invalid sequence data feature in <221> in SEQ ID (2)
W 213	Artificial or Unknown found in <213> in SEQ ID (3)
E 257	Invalid sequence data feature in <221> in SEQ ID (3)
W 213	Artificial or Unknown found in <213> in SEQ ID (4)
E 257	Invalid sequence data feature in <221> in SEQ ID (4)
W 213	Artificial or Unknown found in <213> in SEQ ID (5)
E 257	Invalid sequence data feature in <221> in SEQ ID (5)
W 213	Artificial or Unknown found in <213> in SEQ ID (6)
E 257	Invalid sequence data feature in <221> in SEQ ID (6)

COPY

<210> 1  
<211> 19  
<212> DNA  
<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> forward primer RTGcowM56F2a

<222> (1)..(19)

<223> primer for determination BSE-riskmaterial from cow, sheep and goat

<400> 1

acctgcgacc tggagtcct

19

<210> 2

<211> 15

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> reverse primer RTGcowM56R2a

<222> (1)..(15)

<223> determination of BSE-riskmaterial from cow, goat and sheep

<400> 2

ctcgcgcacg tgccg

15

<210> 3

<211> 17

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

COPY

<220>

<221> TaqManmgb-probe OptiR

<222> (1)..(17)

<223> Determination of BSE riskmaterial from cow, sheep and goat

<400> 3

actcgttcgt gccgcgc

17

<210> 4

<211> 20

<212> DNA  
<213> Artificial  
  
<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> forward primer RTGpigM56F2  
<222> (1)..(20)  
<223> Determination of BSE riskmaterial from pig

<400> 4  
gacctgcgac gtggagtc

20

<210> 5  
<211> 18  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> reverse primer RTGpigM56R2  
<222> (1)..(18)  
<223> Determination of BSE riskmaterial from pig

<400> 5  
tggcgctcct cctgctcc

18

<210> 6  
<211> 17  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> OptiR TaqManmgb probe  
<222> (1)..(17)  
<223> Determination of BSE riskmaterial from pig

<400> 6  
actcggttcgt gccgcgc

17

COPY

# SEQUENCE LISTING

<110> TransMIT Gesellschaft für Technologietransfer mbH

<120> Animal Species-specific and quantitative Detection of Central Nervous  
System Tissue in Meat  
and Meat Products

<130> TM006/Bü/Sc/Ab

<140> DE 103 61 489.3

<141> 2003-12-23

<160> 6

<170> PatentIn version 3.4

<210> 1

<211> 19

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> forward primer RTGcowM56F2a

<222> (1)..(19)

<223> primer for determination BSE-riskmaterial from cow, sheep and  
goat

<400> 1  
acctgcgacc tggagtcct

19

<210> 2

<211> 15

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> reverse primer RTGcowM56R2a

<222> (1)..(15)

<223> determination of BSE-riskmaterial from cow, goat and sheep

<400> 2  
ctcgcgcacg tgccg

15

<210> 3

<211> 17

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

COPY

<220>  
<221> TaqManmgb-probe OptiR  
<222> (1)..(17)  
<223> Determination of BSE riskmaterial from cow, sheep and goat

<400> 3  
actcgttcgt gccgcgc

17

<210> 4  
<211> 20  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> forward primer RTGpigM56F2  
<222> (1)..(20)  
<223> Determination of BSE riskmaterial from pig

<400> 4  
gacctgcgac gtggagtccc

20

<210> 5  
<211> 18  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

COPY

<220>  
<221> reverse primer RTGpigM56R2  
<222> (1)..(18)  
<223> Determination of BSE riskmaterial from pig

<400> 5  
tggcgctcct cctgctcc

18

<210> 6  
<211> 17  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> OptiR TaqManmgb probe  
<222> (1)..(17)  
<223> Determination of BSE riskmaterial from pig

<400> 6  
actcgttcgt gccgcgc

17