

12007-75 Seq List\_English Transl 10-17-08  
SEQUENCE LISTING

<110> TransMIT Gesellschaft fur Technologietransfer mbH

<120> Species-specific and Quantitative Detection of CNS Tissue in Meat  
and Meat Products

<130> TM006/Bü/Sc/Ab

<140> DE 103 61 489.3

<141> 2003-12-23

<160> 6

<170> PatentIn version 3.4

<210> 1

<211> 19

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> forward primer RTGcowM56F2a

<222> (1)..(19)

<223> primer for determination BSE-riskmaterial from cow, sheep and  
goat

<400> 1  
acctgcgacc tggagtcct 19

<210> 2

<211> 15

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> reverse primer RTGcowM56R2a

<222> (1)..(15)

<223> determination of BSE-riskmaterial from cow, goat and sheep

<400> 2  
ctcgcgcac tgccg 15

<210> 3

<211> 17

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

12007-75 Seq List\_English Transl 10-17-08

<221> TaqManmgb-probe OptiR  
 <222> (1)..(17)  
 <223> Determination of BSE riskmaterial from cow, sheep and goat  
  
 <400> 3  
 actcgttcgt gccgcgc 17  
  
 <210> 4  
 <211> 20  
 <212> DNA  
 <213> Artificial  
  
 <220>  
 <223> GFAP-gene Exon 5 and Exon 6  
  
 <220>  
 <221> forward primer RTGpigM56F2  
 <222> (1)..(20)  
 <223> Determination of BSE riskmaterial from pig  
  
 <400> 4  
 gacctgacgtac gtggagtccc 20  
  
 <210> 5  
 <211> 18  
 <212> DNA  
 <213> Artificial  
  
 <220>  
 <223> GFAP-gene Exon 5 and Exon 6  
  
 <220>  
 <221> reverse primer RTGpigM56R2  
 <222> (1)..(18)  
 <223> Determination of BSE riskmaterial from pig  
  
 <400> 5  
 tggcgctcct cctgctcc 18  
  
 <210> 6  
 <211> 17  
 <212> DNA  
 <213> Artificial  
  
 <220>  
 <223> GFAP-gene Exon 5 and Exon 6  
  
 <220>  
 <221> OptiR TaqManmgb probe  
 <222> (1)..(17)  
 <223> Determination of BSE riskmaterial from pig  
  
 <400> 6  
 actcgttcgt gccgcgc 17

# SEQUENCE LISTING

<110> TransMIT Gesellschaft fur Technologietransfer mbH  
 <120> Species-specific and Quantitative Detection of CNS Tissue in Meat  
 and Meat Products  
 <130> TM006/Bü/Sc/Ab  
 <140> DE 103 61 489.3  
 <141> 2003-12-23  
 <160> 6  
 <170> PatentIn version 3.4  
 <210> 1  
 <211> 19  
 <212> DNA  
 <213> Artificial  
 <220>  
 <223> GFAP-gene Exon 5 and Exon 6  
 <220>  
 <221> forward primer RTGcowM56F2a  
 <222> (1)..(19)  
 <223> primer for determination BSE-riskmaterial from cow, sheep and  
 goat  
 <400> 1  
 acctgcgacc tggagtcct 19  
 <210> 2  
 <211> 15  
 <212> DNA  
 <213> Artificial  
 <220>  
 <223> GFAP-gene Exon 5 and Exon 6  
 <220>  
 <221> reverse primer RTGcowM56R2a  
 <222> (1)..(15)  
 <223> determination of BSE-riskmaterial from cow, goat and sheep  
 <400> 2  
 ctcgcgcatc tgccg 15  
 <210> 3  
 <211> 17  
 <212> DNA  
 <213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> TaqManmgb-probe OptiR  
<222> (1)..(17)  
<223> Determination of BSE riskmaterial from cow, sheep and goat

<400> 3  
actcgttcgt gccgcgc

17

<210> 4  
<211> 20  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> forward primer RTGpigM56F2  
<222> (1)..(20)  
<223> Determination of BSE riskmaterial from pig

<400> 4  
gacctgcgac gtggagtccc

20

<210> 5  
<211> 18  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> reverse primer RTGpigM56R2  
<222> (1)..(18)  
<223> Determination of BSE riskmaterial from pig

<400> 5  
tggcgctcct cctgctcc

18

<210> 6  
<211> 17  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> OptiR TaqManmgb probe  
<222> (1)..(17)  
<223> Determination of BSE riskmaterial from pig

<400> 6  
actcgttcgt gccgcgc