

12007-75 Seq List_English Transl 10-17-08
SEQUENCE LISTING

<110> TransMIT Gesellschaft fur Technologietransfer mbH

<120> Species-specific and Quantitative Detection of CNS Tissue in Meat and Meat Products

<130> TM006/Bü/Sc/Ab

<140> DE 103 61 489.3

<141> 2003-12-23

<160> 6

<170> PatentIn version 3.4

<210> 1

<211> 19

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> forward primer RTGcowM56F2a

<222> (1)..(19)

<223> primer for determination BSE-riskmaterial from cow, sheep and goat

<400> 1

acctgcgacc tggagtcct 19

<210> 2

<211> 15

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> reverse primer RTGcowM56R2a

<222> (1)..(15)

<223> determination of BSE-riskmaterial from cow, goat and sheep

<400> 2

ctcgcgcatc tgccg 15

<210> 3

<211> 17

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

12007-75 Seq List_English Transl 10-17-08

<221> TaqManmgb-probe OptiR
 <222> (1)..(17)
 <223> Determination of BSE riskmaterial from cow, sheep and goat

 <400> 3
 actcgttcgt gccgcgc 17

 <210> 4
 <211> 20
 <212> DNA
 <213> Artificial

 <220>
 <223> GFAP-gene Exon 5 and Exon 6

 <220>
 <221> forward primer RTGpigM56F2
 <222> (1)..(20)
 <223> Determination of BSE riskmaterial from pig

 <400> 4
 gacctgcgac gtggagtccc 20

 <210> 5
 <211> 18
 <212> DNA
 <213> Artificial

 <220>
 <223> GFAP-gene Exon 5 and Exon 6

 <220>
 <221> reverse primer RTGpigM56R2
 <222> (1)..(18)
 <223> Determination of BSE riskmaterial from pig

 <400> 5
 tggcgctcct cctgctcc 18

 <210> 6
 <211> 17
 <212> DNA
 <213> Artificial

 <220>
 <223> GFAP-gene Exon 5 and Exon 6

 <220>
 <221> OptiR TaqManmgb probe
 <222> (1)..(17)
 <223> Determination of BSE riskmaterial from pig

 <400> 6
 actcgttcgt gccgcgc 17

SEQUENCE LISTING

<110> TransMIT Gesellschaft für Technologietransfer mbH
 <120> Species-specific and Quantitative Detection of CNS Tissue in Meat
 and Meat Products
 <130> TM006/Bü/Sc/Ab
 <140> DE 103 61 489.3
 <141> 2003-12-23
 <160> 6
 <170> PatentIn version 3.4
 <210> 1
 <211> 19
 <212> DNA
 <213> Artificial
 <220>
 <223> GFAP-gene Exon 5 and Exon 6
 <220>
 <221> forward primer RTGcowM56F2a
 <222> (1)..(19)
 <223> primer for determination BSE-riskmaterial from cow, sheep and
 goat
 <400> 1
 acctgcgacc tggagtcct 19
 <210> 2
 <211> 15
 <212> DNA
 <213> Artificial
 <220>
 <223> GFAP-gene Exon 5 and Exon 6
 <220>
 <221> reverse primer RTGcowM56R2a
 <222> (1)..(15)
 <223> determination of BSE-riskmaterial from cow, goat and sheep
 <400> 2
 ctcgcgcac tgccg 15
 <210> 3
 <211> 17
 <212> DNA
 <213> Artificial

<220>
<223> GFAP-gene Exon 5 and Exon 6

<220>
<221> TaqManmgb-probe OptiR
<222> (1)..(17)
<223> Determination of BSE riskmaterial from cow, sheep and goat

<400> 3
actcgttcgt gccgcgc

17

<210> 4
<211> 20
<212> DNA
<213> Artificial

<220>
<223> GFAP-gene Exon 5 and Exon 6

<220>
<221> forward primer RTGpigM56F2
<222> (1)..(20)
<223> Determination of BSE riskmaterial from pig

<400> 4
gacctgacgt gtggagtccc

20

<210> 5
<211> 18
<212> DNA
<213> Artificial

<220>
<223> GFAP-gene Exon 5 and Exon 6

<220>
<221> reverse primer RTGpigM56R2
<222> (1)..(18)
<223> Determination of BSE riskmaterial from pig

<400> 5
tggcgctcct cctgctcc

18

<210> 6
<211> 17
<212> DNA
<213> Artificial

<220>
<223> GFAP-gene Exon 5 and Exon 6

<220>
<221> OptiR TaqManmgb probe
<222> (1)..(17)
<223> Determination of BSE riskmaterial from pig

<400> 6
actcgttcgt gccgcgc