



# SEQUENCE LISTING

<110> TransMIT Gesellschaft fur Technologietransfer mbH .

<120> Species-specific and Quantitative Detection of CNS Tissue in Meat  
and Meat Products

<130> TM006/Bu/Sc/Ab

<140> DE 103 61 489.3

<141> 2003-12-23

<160> 6

<170> PatentIn version 3.4

<210> 1

<211> 19

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> forward primer RTGcowM56F2a

<222> (1) .. (19)

<223> primer for determination BSE-riskmaterial from cow, sheep and  
goat

<400> 1  
acctgcgacc tggagtcct 19

<210> 2

<211> 15

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> reverse primer RTGcowM56R2a

<222> (1) .. (15)

<223> determination of BSE-riskmaterial from cow, goat and sheep

<400> 2  
ctcgcgcac tgccg 15

<210> 3

<211> 17

<212> DNA

<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> TaqManmgb-probe OptiR  
<222> (1)..(17)  
<223> Determination of BSE riskmaterial from cow, sheep and goat

<400> 3  
actcgttcgt gccgcgc

17

<210> 4  
<211> 20  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> forward primer RTGpigM56F2  
<222> (1)..(20)  
<223> Determination of BSE riskmaterial from pig

<400> 4  
gacctgcgac gtggagtccc

20

<210> 5  
<211> 18  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> reverse primer RTGpigM56R2  
<222> (1)..(18)  
<223> Determination of BSE riskmaterial from pig

<400> 5  
tggcgctcct cctgctcc

18

<210> 6  
<211> 17  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> OptiR TaqManmgb probe  
<222> (1)..(17)  
<223> Determination of BSE riskmaterial from pig

<400> 6  
actcgttcgt gccgcgc



12007-75 Seq List\_English Transl 1-5-09  
SEQUENCE LISTING

<110> TransMIT Gesellschaft fur Technologietransfer mbH  
<120> Species-specific and Quantitative Detection of CNS Tissue in Meat  
and Meat Products  
<130> TM006/Bu/Sc/Ab  
<140> DE 103 61 489.3  
<141> 2003-12-23  
<160> 6  
<170> PatentIn version 3.4  
<210> 1  
<211> 19  
<212> DNA  
<213> Artificial  
<220>  
<223> GFAP-gene Exon 5 and Exon 6  
  
<220>  
<221> forward primer RTGcowM56F2a  
<222> (1)..(19)  
<223> primer for determination BSE-riskmaterial from cow, sheep and  
goat  
  
<400> 1  
acctgcgacc tggagtcct 19  
  
<210> 2  
<211> 15  
<212> DNA  
<213> Artificial  
<220>  
<223> GFAP-gene Exon 5 and Exon 6  
  
<220>  
<221> reverse primer RTGcowM56R2a  
<222> (1)..(15)  
<223> determination of BSE-riskmaterial from cow, goat and sheep  
  
<400> 2  
ctcgcgcatac tgccg 15  
  
<210> 3  
<211> 17  
<212> DNA  
<213> Artificial  
<220>  
<223> GFAP-gene Exon 5 and Exon 6  
  
<220>

12007-75 Seq List\_English Transl 1-5-09

<221> TaqManmgb-probe OptiR  
 <222> (1)..(17)  
 <223> Determination of BSE riskmaterial from cow, sheep and goat  
  
 <400> 3  
 actcgttcgt gccgcgc 17  
  
 <210> 4  
 <211> 20  
 <212> DNA  
 <213> Artificial  
  
 <220>  
 <223> GFAP-gene Exon 5 and Exon 6  
  
 <220>  
 <221> forward primer RTGpigM56F2  
 <222> (1)..(20)  
 <223> Determination of BSE riskmaterial from pig  
  
 <400> 4  
 gacctgacgtac gtggagtcgc 20  
  
 <210> 5  
 <211> 18  
 <212> DNA  
 <213> Artificial  
  
 <220>  
 <223> GFAP-gene Exon 5 and Exon 6  
  
 <220>  
 <221> reverse primer RTGpigM56R2  
 <222> (1)..(18)  
 <223> Determination of BSE riskmaterial from pig  
  
 <400> 5  
 tggcgctcct cctgctcc 18  
  
 <210> 6  
 <211> 17  
 <212> DNA  
 <213> Artificial  
  
 <220>  
 <223> GFAP-gene Exon 5 and Exon 6  
  
 <220>  
 <221> OptiR TaqManmgb probe  
 <222> (1)..(17)  
 <223> Determination of BSE riskmaterial from pig  
  
 <400> 6  
 actcgttcgt gccgcgc 17