



SEQUENCE LISTING

<110> TransMIT Gesellschaft fur Technologietransfer mbH .

<120> Species-specific and Quantitative Detection of CNS Tissue in Meat
and Meat Products

<130> TM006/Bu/Sc/Ab

<140> DE 103 61 489.3

<141> 2003-12-23

<160> 6

<170> PatentIn version 3.4

<210> 1

<211> 19

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> forward primer RTGcowM56F2a

<222> (1)..(19)

<223> primer for determination BSE-riskmaterial from cow, sheep and
goat

<400> 1

acctgcgacc tggagtcct 19

<210> 2

<211> 15

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> reverse primer RTGcowM56R2a

<222> (1)..(15)

<223> determination of BSE-riskmaterial from cow, goat and sheep

<400> 2

ctcgcgcacg tgccg 15

<210> 3

<211> 17

<212> DNA

<213> Artificial

<220>
<223> GFAP-gene Exon 5 and Exon 6

<220>
<221> TaqManmgb-probe OptiR
<222> (1)..(17)
<223> Determination of BSE riskmaterial from cow, sheep and goat

<400> 3
actcgttcgt gccgcgc

17

<210> 4
<211> 20
<212> DNA
<213> Artificial

<220>
<223> GFAP-gene Exon 5 and Exon 6

<220>
<221> forward primer RTGpigM56F2
<222> (1)..(20)
<223> Determination of BSE riskmaterial from pig

<400> 4
gacctgcgac gtggagtcctc

20

<210> 5
<211> 18
<212> DNA
<213> Artificial

<220>
<223> GFAP-gene Exon 5 and Exon 6

<220>
<221> reverse primer RTGpigM56R2
<222> (1)..(18)
<223> Determination of BSE riskmaterial from pig

<400> 5
tggcgctcct cctgctcc

18

<210> 6
<211> 17
<212> DNA
<213> Artificial

<220>
<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> OptiR TaqManmgb probe

<222> (1)..(17)

<223> Determination of BSE riskmaterial from pig

<400> 6

actcgttcgt gccgcgc

17



12007-75 Seq List_English Transl 1-5-09
SEQUENCE LISTING

<110> TransMIT Gesellschaft fur Technologietransfer mbH

<120> species-specific and Quantitative Detection of CNS Tissue in Meat
and Meat Products

<130> TM006/Bu/Sc/Ab

<140> DE 103 61 489.3

<141> 2003-12-23

<160> 6

<170> PatentIn version 3.4

<210> 1

<211> 19

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> forward primer RTGcowM56F2a

<222> (1)..(19)

<223> primer for determination BSE-riskmaterial from cow, sheep and
goat

<400> 1
acctgcgacc tggagtcct 19

<210> 2

<211> 15

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> reverse primer RTGcowM56R2a

<222> (1)..(15)

<223> determination of BSE-riskmaterial from cow, goat and sheep

<400> 2
ctcgcgcac tgccg 15

<210> 3

<211> 17

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

12007-75 Seq List_English Transl 1-5-09

<221> TaqManmgb-probe OptiR
 <222> (1)..(17)
 <223> Determination of BSE riskmaterial from cow, sheep and goat

<400> 3
 actcggttcgt gccgcgc

17

<210> 4
 <211> 20
 <212> DNA
 <213> Artificial

<220>
 <223> GFAP-gene Exon 5 and Exon 6

<220>
 <221> forward primer RTGpigM56F2
 <222> (1)..(20)
 <223> Determination of BSE riskmaterial from pig

<400> 4
 gacctgcgac gtggagtccc

20

<210> 5
 <211> 18
 <212> DNA
 <213> Artificial

<220>
 <223> GFAP-gene Exon 5 and Exon 6

<220>
 <221> reverse primer RTGpigM56R2
 <222> (1)..(18)
 <223> Determination of BSE riskmaterial from pig

<400> 5
 tggcgctcct cctgctcc

18

<210> 6
 <211> 17
 <212> DNA
 <213> Artificial

<220>
 <223> GFAP-gene Exon 5 and Exon 6

<220>
 <221> OptiR TaqManmgb probe
 <222> (1)..(17)
 <223> Determination of BSE riskmaterial from pig

<400> 6
 actcggttcgt gccgcgc

17