

=====

Sequence Listing was accepted.

See attached Validation Report.

If you need help call the Patent Electronic Business Center at (866)
217-9197 (toll free).

Reviewer: Keisha Douglas

Timestamp: [year=2009; month=1; day=22; hr=15; min=44; sec=28; ms=310;]

=====

Application No: 10584270 Version No: 3.0

Input Set:

Output Set:

Started: 2009-01-07 15:30:38.667
Finished: 2009-01-07 15:30:40.008
Elapsed: 0 hr(s) 0 min(s) 1 sec(s) 341 ms
Total Warnings: 6
Total Errors: 6
No. of SeqIDs Defined: 6
Actual SeqID Count: 6

Error code	Error Description
W 213	Artificial or Unknown found in <213> in SEQ ID (1)
E 257	Invalid sequence data feature in <221> in SEQ ID (1)
W 213	Artificial or Unknown found in <213> in SEQ ID (2)
E 257	Invalid sequence data feature in <221> in SEQ ID (2)
W 213	Artificial or Unknown found in <213> in SEQ ID (3)
E 257	Invalid sequence data feature in <221> in SEQ ID (3)
W 213	Artificial or Unknown found in <213> in SEQ ID (4)
E 257	Invalid sequence data feature in <221> in SEQ ID (4)
W 213	Artificial or Unknown found in <213> in SEQ ID (5)
E 257	Invalid sequence data feature in <221> in SEQ ID (5)
W 213	Artificial or Unknown found in <213> in SEQ ID (6)
E 257	Invalid sequence data feature in <221> in SEQ ID (6)

SEQUENCE LISTING

<110> TransMIT Gesellschaft fur Technologietransfer mbH

<120> Species-specific and Quantitative Detection of CNS Tissue in Meat
and Meat Products

<130> TM006/Bu/Sc/Ab

<140> 10584270

<141> 2009-01-07

<150> DE 103 61 489.3

<151> 2003-12-23

<160> 6

<170> PatentIn version 3.4

<210> 1

<211> 19

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> forward primer RTGcowM56F2a

<222> (1)..(19)

<223> primer for determination BSE-riskmaterial from cow, sheep and
goat

<400> 1

acctgcgacc tggagtcct 19

<210> 2

<211> 15

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> reverse primer RTGcowM56R2a

<222> (1)..(15)

<223> determination of BSE-riskmaterial from cow, goat and sheep

<400> 2

ctcgcgcatc tgccg 15

<210> 3

<211> 17
 <212> DNA
 <213> Artificial

 <220>
 <223> GFAP-gene Exon 5 and Exon 6

 <220>
 <221> TaqManmgb-probe OptiR
 <222> (1)..(17)
 <223> Determination of BSE riskmaterial from cow, sheep and goat

 <400> 3
 actcgttcgt gccgcgc 17

<210> 4
 <211> 20
 <212> DNA
 <213> Artificial

 <220>
 <223> GFAP-gene Exon 5 and Exon 6

 <220>
 <221> forward primer RTGpigM56F2
 <222> (1)..(20)
 <223> Determination of BSE riskmaterial from pig

 <400> 4
 gacctgcgac gtggagtccc 20

<210> 5
 <211> 18
 <212> DNA
 <213> Artificial

 <220>
 <223> GFAP-gene Exon 5 and Exon 6

 <220>
 <221> reverse primer RTGpigM56R2
 <222> (1)..(18)
 <223> Determination of BSE riskmaterial from pig

 <400> 5
 tggcgctcct cctgctcc 18

<210> 6
 <211> 17
 <212> DNA
 <213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> OptiR TaqManmgb probe

<222> (1)..(17)

<223> Determination of BSE riskmaterial from pig

<400> 6

actcgttcgt gccgcgc