

=====

Sequence Listing was accepted.

See attached Validation Report.

If you need help call the Patent Electronic Business Center at (866)  
217-9197 (toll free).

Reviewer: Keisha Douglas

Timestamp: [year=2009; month=1; day=22; hr=15; min=44; sec=28; ms=310; ]

=====

Application No: 10584270 Version No: 3.0

**Input Set:****Output Set:**

**Started:** 2009-01-07 15:30:38.667  
**Finished:** 2009-01-07 15:30:40.008  
**Elapsed:** 0 hr(s) 0 min(s) 1 sec(s) 341 ms  
**Total Warnings:** 6  
**Total Errors:** 6  
**No. of SeqIDs Defined:** 6  
**Actual SeqID Count:** 6

Error code	Error Description
W 213	Artificial or Unknown found in <213> in SEQ ID (1)
E 257	Invalid sequence data feature in <221> in SEQ ID (1)
W 213	Artificial or Unknown found in <213> in SEQ ID (2)
E 257	Invalid sequence data feature in <221> in SEQ ID (2)
W 213	Artificial or Unknown found in <213> in SEQ ID (3)
E 257	Invalid sequence data feature in <221> in SEQ ID (3)
W 213	Artificial or Unknown found in <213> in SEQ ID (4)
E 257	Invalid sequence data feature in <221> in SEQ ID (4)
W 213	Artificial or Unknown found in <213> in SEQ ID (5)
E 257	Invalid sequence data feature in <221> in SEQ ID (5)
W 213	Artificial or Unknown found in <213> in SEQ ID (6)
E 257	Invalid sequence data feature in <221> in SEQ ID (6)

# SEQUENCE LISTING

<110> TransMIT Gesellschaft fur Technologietransfer mbH

<120> Species-specific and Quantitative Detection of CNS Tissue in Meat  
and Meat Products

<130> TM006/Bu/Sc/Ab

<140> 10584270

<141> 2009-01-07

<150> DE 103 61 489.3

<151> 2003-12-23

<160> 6

<170> PatentIn version 3.4

<210> 1

<211> 19

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> forward primer RTGcowM56F2a

<222> (1)..(19)

<223> primer for determination BSE-riskmaterial from cow, sheep and  
goat

<400> 1

acctgcgacc tggagtcct 19

<210> 2

<211> 15

<212> DNA

<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> reverse primer RTGcowM56R2a

<222> (1)..(15)

<223> determination of BSE-riskmaterial from cow, goat and sheep

<400> 2

ctcgcgcac tgcgcg 15

<210> 3

<211> 17  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> TaqManmgb-probe OptiR  
<222> (1)..(17)  
<223> Determination of BSE riskmaterial from cow, sheep and goat

<400> 3  
actcgttcgt gccgcgc 17

<210> 4  
<211> 20  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> forward primer RTGpigM56F2  
<222> (1)..(20)  
<223> Determination of BSE riskmaterial from pig

<400> 4  
gacctgcgac gtggagtcce 20

<210> 5  
<211> 18  
<212> DNA  
<213> Artificial

<220>  
<223> GFAP-gene Exon 5 and Exon 6

<220>  
<221> reverse primer RTGpigM56R2  
<222> (1)..(18)  
<223> Determination of BSE riskmaterial from pig

<400> 5  
tggcgctcct cctgctcc 18

<210> 6  
<211> 17  
<212> DNA  
<213> Artificial

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> OptiR TaqManmgb probe

<222> (1)..(17)

<223> Determination of BSE riskmaterial from pig

<400> 6

actcggttcgt gccgcgc